

Waterline Maintenance Tablets

BluTab® Waterline Maintenance Tablets are formulated specifically for dental unit waterlines. BluTab® is designed to be continuously present in your waterlines and helps keep lines clean.

Directions for Use:

- 1) Place tablet or tablets (see table below) in empty dental unit water bottle every time you refill the water. BluTab® is intended for use with potable water (including tap or distilled).
- 2) Fill bottle with desired volume of water and connect to dental unit. Tablet dissolves in a few minutes.

Water Volume	Tablets
2 liters (68 fl.oz)	1
4 liters (135 fl. oz.)	2

- 3) BluTab® may be left in the waterlines overnight and for extended periods of non-use. BluTab® has been shown to be effective for up to 28 days. Use water solution until it is expended; or up to 28 days. It is not necessary to purge lines at the end of each day.
- 4) BluTab® is not intended for use as a biofilm removal agent. Before starting a treatment program with BluTab®, flush lines with an antimicrobial product registered with the U.S. EPA which is appropriate for use in dental unit water lines. Repeat flushing protocol every two months. Consult with the manufacturer of your dental unit or water delivery system to determine the best method to flush your water system. For a list of recommended products, call 1-800-819-3336.

Clinical monitoring of water quality can ensure that procedures are performed correctly. Test water quality monthly to establish effectiveness of your maintenance program. If water quality results approach or exceed the CDC recommended bacteria counts, repeat flushing protocol.

Precautionary Statement:

Caution: Eye: Tablet may cause irritation to the eyes (diluted product is non-irritating). Avoid contact with skin, eyes or clothing.
Skin: Prolonged contact may cause irritation and redness. Wash hands with soap and water after handling.

First Aid - Eyes: Immediately flush eyes with water for at least 10 minutes. If irritation persists get medical attention.

Skin: For skin contact flush with plenty of water. If irritation persists get medical attention. Take off any contaminated clothing; wash before reuse.



109 Inverness Drive East, Unit F
 Englewood, CO 80112-5105
 Toll Free: 1-800-819-3336

Waterline Maintenance Tablets

BluTab® Waterline Maintenance Tablets are formulated specifically for dental unit waterlines. BluTab® is designed to be continuously present in your waterlines and helps keep lines clean.

Directions for Use:

- 1) Place tablet or tablets (see table below) in empty dental unit water bottle every time you refill the water. BluTab® is intended for use with potable water (including tap or distilled).
- 2) Fill bottle with desired volume of water and connect to dental unit. Tablet dissolves in a few minutes.

Water Volume	Tablets
700-750 mls (24-25 fl.oz.)	1
1500 mls (51 fl.oz.)	2
2000 mls (68 fl.oz.)	3

- 3) BluTab® may be left in the waterlines overnight and for extended periods of non-use. BluTab® has been shown to be effective for up to 28 days. Use water solution until it is expended; or up to 28 days. It is not necessary to purge lines at the end of each day.
- 4) BluTab® is not intended for use as a biofilm removal agent. Before starting a treatment program with BluTab®, flush lines with an antimicrobial product registered with the U.S. EPA which is appropriate for use in dental unit water lines. Repeat flushing protocol every two months. Consult with the manufacturer of your dental unit or water delivery system to determine the best method to flush your water system. For a list of recommended products, call 1-800-819-3336.

Clinical monitoring of water quality can ensure that procedures are performed correctly. Test water quality monthly to establish effectiveness of your maintenance program. If water quality results approach or exceed the CDC recommended bacteria counts, repeat flushing protocol.

Precautionary Statement:

Caution: Eye: Tablet may cause irritation to the eyes (diluted product is non-irritating). Avoid contact with skin, eyes or clothing.
Skin: Prolonged contact may cause irritation and redness. Wash hands with soap and water after handling.

First Aid - Eyes: Immediately flush eyes with water for at least 10 minutes. If irritation persists get medical attention.

Skin: For skin contact flush with plenty of water. If irritation persists get medical attention. Take off any contaminated clothing; wash before reuse.



109 Inverness Drive East, Unit F
 Englewood, CO 80112-5105
 Toll Free: 1-800-819-3336