

ConFirm®

Complete In-Office Biological Monitoring System



Results within 46 hours!

Fast Results

Now you can get results in 24 hours with ConFirm's In-Office Biological Monitoring System.

Easy to Use

Simply process a vial in a normal sterilizer cycle. Activate the vial, place it in the incubator, and read and record results.

Easy to Read Results

If any bacterial spores survive the process, the vial media will turn yellow.

Economical

ConFirm In-Office System is one of the most cost effective methods of biological monitoring.

Complete

Everything you need for effective biological monitoring of steam or ethylene oxide sterilizers.

Universal

ConFirm Biological Indicators can be processed in any appropriate incubator (55-60°C for steam, 35-39° C for gas), including the current 3M Attest® Incubator. Attest is a registered trademark of 3M Corp.

Order Information

Description	Order Number
Steam Biological Indicators 25/box	CSB125
Incubator 14 well, dual temperature	CIN14
Steam System 25 indicators, incubator, crusher and record keeper	CBMS10
Ethylene Oxide Biological Indicators 25/box	CSCG25
Ethylene Oxide System 25 indicators, incubator, crusher and record keeper	CSCGS



109 Inverness Drive East, Unit F
Englewood, Colorado 80112-5105
Toll Free: 1-800-819-3336
www.confirmmonitoring.com
L-CI001