

# ConFirm®

## Complete In-Office Biological Monitoring System



**Results within 48 hours!**

### **Fast Results**

Now you can get results in 24 to 48 hours with ConFirm's In-Office Biological Monitoring System.

### **Easy to Use**

Simply process a vial in a normal sterilizer cycle. Activate the vial, place it in the incubator, and read and record results.

### **Easy to Read Results**

If any bacterial spores survive the process, the vial media will turn yellow.

### **Economical**

ConFirm In-Office System is one of the most cost effective methods of biological monitoring.

### **Complete**

Everything you need for effective biological monitoring of steam or ethylene oxide sterilizers.

### **Universal**

ConFirm Biological Indicators can be processed in any appropriate incubator (55-60°C for steam, 35-39° C for gas), including the current 3M Attest® Incubator. Attest is a registered trademark of 3M Corp.

#### Order Information

Description	Order Number
Steam Biological Indicators 25/box	CSB125
Incubator 14 well, dual temperature	CIN14
Steam System 25 indicators, incubator, crusher and record keeper	CBMS10
Ethylene Oxide Biological Indicators 25/box	CSCG25
Ethylene Oxide System 25 indicators, incubator, crusher and record keeper	CSCGS

 **ConFirm®**  
Monitoring Systems, Inc.  
14550 E. Easter Avenue  
Suite 600  
Englewood, CO 80112-4263  
Toll Free: 1-800-819-3336  
[www.confirmmonitoring.com](http://www.confirmmonitoring.com)